



MediLumine Inc.
5795 Av. De Gaspé, Suite #112
Montreal, Quebec, Canada
H2S 2X3

Product Information Sheet


Description

Product: Bradykinin B1 Receptor Agonist	Quantity: 1 mg
Sequence: NH₂-Sar-Lys-Arg-Pro-Pro-Gly-Phe-Ser-Pro-(D-Phe)-COOH	Production date: August 2015
Application: Permeabilization of BBTB barrier	Storage condition: -20 degrees Celcius and dryness (stable for 1 year)
Cas#: 1158908-54-3	Long term storage: -70 degrees Celcius in dark and dryness

Final product is in lyophilized powder form.

Positive

MS Spectrum



Sample Information Item

Date and Time :2015-7-17 10:15:21
Usr :
Sample :MEDL003-1 Sar-FK-10
Inj. Volume 1
MW :1103.30

Probe: ESI
Nebulizer Gas Flow: 1.5L/min
CDL: -20.0v
DL Temp: 250°C
Block Temp: 300°C

Probe bias: +4.5kv
Detector: 1.5kv
T. Flow: 0.2ml/min
B. conc: 50%H2O/50%ACN